Instructions for Use

REV.08

Effective Date: May. 30, 2019

IFU for Tytrak PTCA balloon catheter	1
Mode d'emploi pour Cathéter à ballonnet PTCA Tytrak	4
Istruzioni per l'uso del catetere a palloncino Tytrak per PTCA	7
Instrucciones De Uso Del Catéter Con Balón Para ACTP Tytrak	10
Gebrauchsanweisung für den Tytrak-PTCA-Ballonkatheter	13
Gebbruiksaanwijzing voor Tytrak PTCA-ballonkatheter	16
Instruções de Utilização Do Tytrak – Cateter Balão PTCA	19
Tytrak PTCA 氣球導管使用說明	22

0.009.0083

SINOMED
Innovation for health

Single Use Only

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Esclusivamente monouso

Dispositivo de un solo uso

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Fabrikant

Fabricante

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Authorised representative in the EU

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IFU for Tytrak PTCA Balloon Catheter

DEVICE NAME

- Tytrak PTCA Balloon Catheter (Trade Mark: Tytrak)

PACKING LIST

- A set of Tytrak PTCA Balloon Catheter
- A copy of IFU

WARNING

- Before use, please read the IFU carefully and pay more attention to all the items in Warning and Precautions to keep from accidental events.
- Selection of patients is necessary since the use of this device carries the risks associated with PTCA, including acute occlusion of the artery; thrombosis (acute or subacute); dissection; hemorrhage or bleeding events.
- Tytrak PTCA Balloon Catheter is sterilized with ETO and non-pyrogenic. For single use only. Any reuse, reprocessing or resterilization may result in the structural deformation, failure, contamination of the device which create a potential risk of, cross-infection. What mentioned above can lead injury, disease or death.
- Only the qualified physicians who have received professional training are permitted to perform the surgery of the balloon catheter.
- Before use, please check the package carefully. Do not use if the packaging box or the inner packaging bag is opened or damaged.
- Before use, NOTE the "Use By" date on the product label carefully, do not use the product over the expiration date.
- Do not exceed rated burst pressure (RBP) indicated on the product label.


1. DEVICE DESCRIPTIONS

Tytrak PTCA Balloon Catheter is a rapid exchange catheter. The Balloon Diameter ranges from 1.25mm to 4.00mm and the Balloon Length ranges from 6mm to 35mm.

Tytrak balloon catheter is used for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis to improve myocardial perfusion.

2. DEVICE STRUCTURE

Tytrak PTCA Balloon Catheter is composed of Polypropylene, PTFE, Pebax, Polyolefin, Stainless steel, etc.

- | | | |
|-------------------|------------------------|-------------------------------------|
| 1) Balloon TIP | 5) Marker | 9) Luer Lock |
| 2) Balloon Sheath | 6) Outer Tubing | Note 1: D1 is balloon diameter |
| 3) Balloon | 7) Hypotube | Note 2: L2 is balloon length |
| 4) Inner Tubing | 8) Luer Lock Protector | Note 3: L is delivery system length |

There are two delivery system shaft markers (93 cm and 103cm from the distal tip) , which indicate the relative position of the delivery system to the end of a brachial or femoral

guiding catheter, and two radiopaque markers, located underneath the balloon, which fluoroscopically mark the working length of the balloon. And near the tip there is a notch for the guide wire which is 24 cm from the distal tip. The guiding catheter with inner diameter 1.42 mm (5Fr/0.056") and the guide wire with diameter 0.36 mm (0.014") is suitable for Tytrak PTCA Balloon Catheter.

3. DEVICE SPECIFICATIONS

SPECIFICATIONS OF BALLOON CATHETER

Diameter (mm)	Length (mm)							
	6	10	12	15	20	25	30	35
1.25		●		●	●			
1.50		●		●	●			
2.00		●		●	●	●		
2.25		●		●	●	●	●	●
2.50	●	●	●	●	●	●	●	●
2.75	●	●	●	●	●	●	●	●
3.00	●	●	●	●	●	●	●	●
3.25	●	●	●	●	●	●	●	●
3.50	●	●	●	●	●	●	●	●
4.00	●	●	●	●	●	●	●	●

NOMINAL PRESSURE AND RATED BUST PRESSURE OF BALLOON CATHETER

Diameter (mm)	Inflation Pressure (atm)						
	4	6	8	10	12	14	16
1.25	1.18	1.25	1.28	1.32	1.36	1.40	1.44
1.50	1.42	1.50	1.53	1.56	1.59	1.62	1.65
2.00	1.87	2.00	2.06	2.12	2.17	2.21	2.26
2.25	2.19	2.25	2.29	2.35	2.39	2.43	2.47
2.50	2.32	2.50	2.59	2.66	2.73	2.79	2.84
2.75	2.48	2.75	2.81	2.88	2.94	3.00	3.06
3.00	2.84	3.00	3.08	3.17	3.24	3.31	3.38
3.25	3.18	3.25	3.32	3.39	3.45	3.52	3.57
3.50	3.37	3.50	3.58	3.65	3.73	3.81	3.89
4.00	3.85	4.00	4.09	4.19	4.29	4.38	4.47

■ NP ■ RBP

Notes: The rated burst pressure is based on the results of in vitro testing. At least 99.9

percent of balloons (with a 95 percent confidence) will not burst at or below this pressure. Using a pressure monitoring device is recommended to prevent over pressurization. The size of balloon expansion should be confirmed by radiography.

4. INDICATIONS

Tytrak PTCA Balloon Catheter is used for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis to improve myocardial perfusion.

5. CONTRAINDICATIONS

- Patients with all the contraindications of angioplasty;
- Patients with a hypersensitivity to contrast medium;
- Patients with a hypersensitivity to Pebax or stainless steel 304;
- Patients who have dysfunction of blood coagulation or bleeding tendency;
- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated;
- Patients with serious calcified lesions;
- Other patients judged to have a lesion that prevents complete inflation of an angioplasty balloon.

6. COMPLICATIONS

ADVERSE EVENTS, WHICH MAY BE ASSOCIATED WITH THE INTERVENTION OF A BALLOON CATHETER IN CORONARY ARTERIES, BUT NOT LIMITED TO THE FOLLOWING:

- Allergic reaction;
- Arrhythmias;
- Aneurysm;
- Hemorrhage;
- Hypotension/Hypertension;
- Death;
- Dissection;
- Vessel spasm;
- Infection and pain at the insertion site.

PRECAUTIONS

A. Individualization of Treatment

The risks and benefits should be considered for each patient before use of the device.

B. Use in Special Populations

① Pregnancy

There are no adequate and well-controlled studies in pregnant women of Tytrak Balloon Catheter. Effective contraception should be initiated before use of a Tytrak Balloon Catheter. The Tytrak Balloon Catheter should be used during pregnancy only if the potential benefit exceeds the potential risk to the embryo or fetus.

② Lactation

A decision should be made whether to discontinue nursing or to use of the balloon, taking into account the importance of the intervention to the mother.

③ Pediatric use

The safety and efficacy of the Tytrak Balloon Catheter in pediatric patients below the age of 18 have not been testified.

④ Geriatric use

There are no enough studies can find clinically significant differences in outcomes between patients under 65 and over 65 years of age treated with Tytrak PTCA Balloon Catheter.

7. OPERATOR'S MANUAL

ACCESSORIES REQUIRED

- Appropriate Guiding Catheter;
- Contrast Catheter;
- Percutaneous puncturing needle and short guide wire;
- Syringe (10~ 20CC);
- Heparinized Normal Saline and Contrast Medium;
- Inflation Device;
- Three-way stopcock;
- 0.014" (0.36mm) (or less diameter) Guide wire and Hole punch holder of guide wire.

7.1 INSPECTION PRIOR TO USE

DO NOT USE THE DEVICE UNDER THE CONDITIONS NOTED BELOW:

- Do not use if the expiration date has over passed;
- Check the package. Do not use if the packaging box or the inner packaging bag is opened or damaged.

PRECAUTIONS

- Before use, please read this IFU carefully and pay more attention to all the items in Warning and Precautions to keep from accidental events ;
- Only the qualified physicians who have received professional training are permitted to perform the surgery of the balloon catheter implantation;
- For single use only. Do not resterilize or reuse.

7.2 UNPACKAGING

- Tear open the inner pouch to remove the product carefully;
- The procedure of removing avoids bends, kinks, and other damages;
- Carefully check the product to make sure there are no bends, kinks, and other damages ;
- Remove the protecting sheath and protecting wire.

PRECAUTIONS

- Only the appropriate balloon inflation method should be use. Do not use air or any gaseous medium to inflate the balloon in order to avoid the expansion and difficulty in delivery of the balloon;
- The Balloon diameter and length marked on the label means that the diameter and length of the balloon inflated under the rated pressure.

7.3 PREPARATION

- Prepare a syringe filled with heparinized normal saline, attach the syringe to the flushing tool and insert the flushing tool into the distal end of the catheter, inject heparinized normal saline into the lumen. Repeat flushing procedure for several times. Heparinized normal saline should exit from the notch.

PRECAUTIONS

- Do not inflate or deflate the balloon during the preparation;
- Submerge the balloon in sterile heparinized normal saline during balloon preparation to activate the coating.

7.4 DELIVERY PREPARATION

- Prepare the inflation device or syringe with suitable diluted contrast medium;
- Attach the inflation device or syringe to the stopcock; attach to the balloon inflation port hub;
- Keep the pump under medium pressure.

PRECAUTIONS

- Balloon Catheter shall be performed as a single unit.

7.5 DELIVERY PROCEDURE

1. Prepare the vascular access site according to standard practice;
2. Backload the balloon catheter onto the proximal portion of the guidewire while maintaining the guidewire position across the target lesion;
3. Advance the balloon catheter over the guidewire to the target lesion. Use the radiopaque balloon markers to position the balloon across the lesion; Perform angiography to confirm the position of the balloon.

PRECAUTIONS

- Do not expand the balloon if it is not properly positioned in the vessel;
- Should unusual resistance be felt at any time during either lesion access or removal of the balloon catheter before balloon inflation, the entire balloon catheter should be removed as a single unit.

7.6 INFLATION PROCEDURE

1. Before inflation, reconfirm the correct position of the balloon relative to the target lesion via the radiopaque balloon markers;
2. Attach the inflation device (only partially filled with contrast media) to a three-way stopcock and apply negative pressure to purge the balloon of air;
3. Turn the stopcock on the catheter to the off position and purge the inflation device of air. Close the side port of the stopcock;
4. Under fluoroscopic visualization, inflate the balloon to at least the nominal pressure, but do not exceed the labeled rated burst pressure;
5. Deflate the balloon by pulling a vacuum with the inflation device. Make certain that the balloon is fully deflated before attempting to move the catheter.

PRECAUTIONS

- Do not expand the balloon if it is not properly positioned in the vessel;
- Place of the balloon has the potential to compromise side branch patency;
- Balloon pressure should be monitored during inflation. Do not exceed rated burst pressure as indicated on the product label, or they may result in a ruptured balloon with possible intimal damage and dissection;
- Should unusual resistance be felt at any time during either lesion access or removal of the balloon catheter, the entire system should be removed as a single unit.

7.7 REMOVAL PROCEDURE

1. Ensure that the balloon is fully deflated;
2. While maintaining the guidewire position and negative pressure on the inflation device, withdraw the balloon catheter;
3. Repeat angiography to assess the treated area.

PRECAUTIONS

- Should unusual resistance be felt at any time during either lesion access or removal of the balloon catheter, the entire system should be removed as a single unit;
- Tighten the rotating haemostatic valve to secure the balloon catheter to the guiding catheter; then remove the guiding catheter and balloon catheter as a single unit.

8. STERILE

This device is sterilized with ETO and is non-pyrogenic.

9. STORE

AVOID the stress and impact. Keep from direct sunlight or corrosive gas. Store the products in a clean, breezy place, the relative humidity shall be lower than 80%.

Store this device at room temperature (20±10°C).

10. USE BY

The product shelf life is valid for two years from the date of sterilization.

11. WARRANTY AND LIMITATIONS

Sino Medical Sciences Technology Inc. has committed itself to the manufacture of the Tytrak PTCA Balloon Catheter. Sino Medical Sciences Technology Inc. warrants that all the products shall be free of defects in materials and manufacturing.

12. REUSE PRECAUTION STATEMENT

Do not use if sterile barrier is damaged.

For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

13. DECLARATION

Sino Medical Sciences Technology Inc. will not to be liable for any accident and damage result from reuse of the product, incorrect size selection of the balloon, misoperation or any other mistakes caused by human-factor(directly or indirectly).

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SINOMED B.V.