HT SUPREME™

Drug Coated Coronary Stent System

Instruction for Use

REV: 10

Effective Date: 04 March, 2021

English	HT Supreme™ Drug Coated Coronary Stent System	Instruction for Use	1	en
Français	Système d'endoprothèse coronaire à revêtement médicamenteux HT Supreme™	Mode d'emploi	6	fr
Italiano	Sistema stent coronarico rivestito di farmaco HT Supreme™	Istruzioni per l'uso	12	it
Español	Sistema de stent coronario recubierto de fármaco HT Supreme™	Instrucciones de uso	18	es
Deutsch	HT Supreme™ medikamentenbeschichtetes Koronarstentsystem	Verwendungshinweise (IFU)	24	de
Nederlands	HT Supreme™ coronair stentsysteem met geneesmiddelcoating	Gebruiksaanwijzing	31	nl
Português	Sistema de Stent Coronário Revestido por Fármaco HT Supreme™	Instruções de Utilização	37	pt
繁體中文	HT Supreme™ 藥物塗層冠狀動脈支架系統	使用說明	43	cn
Bahasa Indonesia	HT Supreme™ Sistem Stent Koroner Berlapis Obat	Petunjuk Penggunaan	47	id

Do not re-use
Ne pas réutiliser
Non riutilizzare
No reutilizar
Ausschließlich zur einmaligen Verwendung
Niet hergebruiken
Não reutilizar
請勿重複使用
Jangan menggunakan ulang

Use-by date
Date limite de consommation
Data di scadenza
Fecha de caducidad
Ablaufdatum
Te gebruiken tot
Não utilizar após (data)
使用期限
Tanggal kedaluwarsa

Batch code
Code du lot
Codice di lotto
Código de lote
Chargennummer
Batchcode
Código do lote
批次代碼
Kode batch

Content one
Contenu Un
Contenuto uno
Contenido único
Inhalt: ein Stück
Inhoud één
Contem um
內含一件物品
Isi satu

Serial number
Numéro de série
Numero di serie
Número de serie
Seriennummer
Seriennummer
Número de série
序號
Nomor seri

Consult instructions for use
Consulter le mode d'emploi
Consultare le istruzioni per l'uso
Consulte las instrucciones de uso
Verwendungshinweise beachten
Raadpleeg de gebruiksaanwijzing
Consultar as instruções de utilização
請參考使用說明
Lihat petunjuk penggunaan

Do not use if package is damaged
Ne pas utiliser si l'emballage est endommagé
Non utilizzzare se la confezione appare danneggiata.
No utilizar si el envase se encuentra dañado
Nicht verwenden, falls Verpackung beschädigt ist
Niet gebruiken als de verpakking beschadigd is
Não utilizar se a embalagem estiver danificada
如包裝已損壞，請勿使用
Jangan digunakan jika kemasan rusak


Sterilized using irradiation
Stérilisé par irradiation
Sterilizzato con irradiazione
Esterilizado mediante radiación
Durch Bestrahlung sterilisiert
Gesteriliseerd met straling
Esterilizado por irradiação
輻射滅菌
Distrikan menggunakan iradiasi

MR conditional
Compatible avec l'IRM
MR condizionale
Compatible con RM
Bedingt MR-kompatibel
MRI-veilig onder bepaalde voorwaarden
Condicional para RM
與核磁共振成像相容
MR conditional

Non-pyrogenic
Apyrogène
Apirogeno
Apirógeno
Nicht pyrogen
Niet-pyrogeen
Não pirogénico
無熱原
Non-pirogenik

Temperature limit
Limite de température
Limite di temperatura
Limite de temperatura
Temperaturgrenzwert
Temperatuurlimiet
Limite de temperatura
溫度限制
Batas suhu

Manufacturer
Fabricant
Produttore
Fabricante
Hersteller
Fabrikant
Fabricante
製造商
Produsen

Authorized representative in the European Community
Représentant agréé pour la Communauté européenne
Rappresentante autorizzato nella Comunità Europea
Representante autorizado en la Comunidad Europea
Bevollmächtigter Vertreter in der Europäischen Gemeinschaft
Geautoriseerde vertegenwoordiger in de Europese Unie
Representante autorizado na União Europeia
歐盟授權代表
Perwakilan resmi di Komunitas Eropa

Keep Dry
Garder au sec
Mantenere asciutto
Mantener seco
Vor Nässe schützen
Droog houden
Manter seco
保持乾燥
Jaga Agar Tetap Kering

Protect from heat and radioactive sources
Protéger des sources de chaleur et de radioactivité
Proteggere da calore e sorgenti radioattive
Proteger contra las fuentes de calor y de radiación
Vor Wärme- und radioaktiven Quellen schützen
Beschermen tegen hitte en radioactieve bronnen
Proteger do calor e de fontes radioativas
防止熱源和輻射源
Lindungi dari panas dan sumber radioaktif

Do not re-sterilize
Ne pas restériliser
Non risterilizzare
No reesterilizar
Nicht resterilisieren
Niet opnieuw steriliseren
Não reesterilizar
不能二次滅菌
Jangan resterilize

en

fr

it

es

de

nl

pt

cn

id

1. DEVICE DESCRIPTION

The HT Supreme™ Drug Coated Coronary Stent System (HT Supreme™ stent system) is a device/drug combination product consisting of a drug-coated balloon expandable cobalt chromium (CoCr) alloy coronary stent, and a delivery system. The stent is coated with a non-erodible base polymer layer and a drug eluting biodegradable top polymer layer. The drug eluting biodegradable polymer layer consists of a blend of sirolimus and biodegradable polymer. The stent is delivered using a rapid exchange balloon expandable delivery system. The main characteristics of the HT Supreme™ stent system are provided in Table 1.1.

Table 1.1: HT Supreme™ Stent System Product Component Description

Parameter	Characteristics
Stent Characteristics	
Available Stent Lengths (mm)	10, 15, 20, 25, 30, 35
Available Stent Diameters (mm)	2.25, 2.50, 2.75, 3.00, 3.25, 3.50, 4.00
Stent Material	Medical grade L-605 cobalt chromium (CoCr) alloy
Base Polymer Layer	A non-erodible layer of poly n-butyl methacrylate (PBMA) covalently bonded to the stent surface.
Drug Component	A conformal coating of a biodegradable polymer loaded with 1.2 µg/mm ² of sirolimus with a maximum nominal drug content of 285 µg on the large stent 4.00 x 35 mm).
Delivery System Characteristics	
Working Length (cm)	145
Design	<ul style="list-style-type: none"> Single access port to inflation lumen. Guide wire exit notch is located 23.5 cm from tip. Designed for guide wires ≤ 0.014" (0.36 mm). Two proximal delivery system shaft markers (87 cm and 97 cm proximal to the distal tip) indicate the relative position of the delivery system to the end of the brachial or femoral artery guiding catheter.
Balloon Component	A semi-compliant balloon with two radiopaque markers located on the catheter shaft to indicate balloon positioning and expanded stent length.
Balloon Inflation Pressure	<ul style="list-style-type: none"> Nominal inflation pressure: 10 atm (1013 kPa) for 2.25, 2.50, 2.75, 3.00, 3.50 diameters and 12 atm (1216 kPa) for 4.00 diameter Rated Burst Pressure (RBP): 18 atm (1824 kPa) for 2.25, 2.50, 2.75, 3.00, 3.25 and 3.50 diameters and 16 atm (1621 kPa) for 4.00 diameter
Guiding Catheter Inner Diameter	5 F, ≥0.056" I.D. (1.42 mm)
Catheter Shaft Outer Diameter (nominal)	<ul style="list-style-type: none"> Proximal: 2.0 F (0.68 mm) Distal: 2.9 F (0.95 mm)
System Crossing Profile	<ul style="list-style-type: none"> Distal: 3.6F (≤1.20mm) for 2.25, 2.50, 2.75, 3.00mm diameter 3.9F (≤1.30mm) for 3.25, 3.50, 4.00 diameter Proximal: 2.9 F (≤0.95mm) for all sizes

1.1 Device Component Description

The HT Supreme™ stent is composed of L605 cobalt chromium (CoCr) alloy and is available in multiple diameter and length combinations as described in Table 1.1. The stent is coated with a thin layer of non-erodible polymer, poly n-butyl methacrylate (PBMA), that is covalently bound to the stent surface. A top coat consists of sirolimus embedded in poly(lactic co-glycolic acid (PLGA) which is biodegradable. The HT Supreme™ stent is delivered using a rapid exchange balloon expandable delivery system with a working length of 145 cm, compatible with guide wires ≤ 0.014" and guide catheter inner diameters ≥ 5 F (0.056"). Radiopaque markers located on the catheter shaft indicate the working length of the balloon.

1.2 Drug and Polymer Component Descriptions

The HT Supreme™ stent is coated with the drug, sirolimus (active ingredient), mixed in a biodegradable polymer (inactive ingredient).

A base polymer layer (inactive ingredient) is designed to secure adhesion of the drug/biodegradable polymer top coating to the stent.

1.2.1 Sirolimus (Active Ingredient)

The active pharmaceutical ingredient in the HT Supreme™ stent is sirolimus, also known as Rapamycin. Sirolimus is a macrocyclic lactone produced by *Streptomyces hygroscopicus*. The chemical name of sirolimus is (3S, 6R, 7E, 9R, 10R, 12R, 14S, 15E, 17E, 19E, 21S, 23S, 26R, 27R, 34aS)-9, 10, 12, 13, 14, 21, 22, 23, 24, 25, 26, 27, 32, 33, 34, 34a -hexadecahydro-9, 27 - dihydroxy-3- [(1R)-10,21-dimethoxy-4-hydroxy-3-methoxycyclohexyl] -1- methylethyl]-2-pyrrolo [2, 1-c] [1,4] Oxazacycloheptriacontane-1, 5, 11, 28, 29 (4H, 6H, 31H)-pentone. Its molecular formula is C₅₁H₇₉NO₁₃ and its molecular weight is 914.2 (g/mol). The chemical structure is shown in Figure 1.1.


Figure 1.1 Sirolimus Chemical Structure

1.2.2 Base Polymer Layer and Biodegradable Polymer (Inactive Ingredients)

The HT Supreme™ stent contains inactive ingredients including non-erodible poly butyl-methacrylate (PBMA), a base polymer layer, and biodegradable poly (DL-lactic-co-glycolide) (PLGA), the drug matrix layer that contains sirolimus. PBMA is a homopolymer electro-grafted onto the stent surface. The PLGA drug matrix copolymer is mixed with sirolimus and applied to the PBMA coated stent. PLGA is a biodegradable linear aliphatic polymeric polymer. The chemical structures are shown in Figure 1.2.


Figure 1.2: Polymer Coatings Chemical Structures

1.2.3 Product Matrix and Sirolimus Content

Table 1.2 HT Supreme™ Stent System Product Matrix and Sirolimus Content

Catalog Number	HT Supreme™ Component			Delivery Catheter
	Expanded Stent Diameter (mm)	Stent Length (mm)	Nominal Sirolimus Content (µg)	Balloon Size (mm × mm)
8.007.1805	2.25	10	59	2.25 × 10
8.007.1806	2.25	15	89	2.25 × 15
8.007.1807	2.25	20	120	2.25 × 20
8.007.1808	2.25	25	150	2.25 × 25
8.007.1809	2.25	30	180	2.25 × 30
8.007.1810	2.25	35	210	2.25 × 35
8.007.1811	2.50	10	59	2.50 × 10
8.007.1812	2.50	15	89	2.50 × 15
8.007.1813	2.50	20	120	2.50 × 20
8.007.1814	2.50	25	150	2.50 × 25
8.007.1815	2.50	30	180	2.50 × 30
8.007.1816	2.50	35	210	2.50 × 35
8.007.1817	2.75	10	59	2.75 × 10
8.007.1818	2.75	15	89	2.75 × 15
8.007.1819	2.75	20	120	2.75 × 20
8.007.1820	2.75	25	150	2.75 × 25
8.007.1821	2.75	30	180	2.75 × 30
8.007.1822	2.75	35	210	2.75 × 35
8.007.1823	3.00	10	59	3.00 × 10
8.007.1824	3.00	15	89	3.00 × 15
8.007.1825	3.00	20	120	3.00 × 20
8.007.1826	3.00	25	150	3.00 × 25

Catalog Number	HT Supreme™ Component			Delivery Catheter
	Expanded Stent Diameter (mm)	Stent Length (mm)	Nominal Sirolimus Content (µg)	Balloon Size (mm × mm)
8.007.1827	3.00	30	180	3.00 × 30
8.007.1828	3.00	35	210	3.00 × 35
8.007.1829	3.25	10	80	3.25 × 10
8.007.1830	3.25	15	121	3.25 × 15
8.007.1831	3.25	20	162	3.25 × 20
8.007.1832	3.25	25	203	3.25 × 25
8.007.1833	3.25	30	244	3.25 × 30
8.007.1834	3.25	35	285	3.25 × 35
8.007.1835	3.50	10	80	3.50 × 10
8.007.1836	3.50	15	121	3.50 × 15
8.007.1837	3.50	20	162	3.50 × 20
8.007.1838	3.50	25	203	3.50 × 25
8.007.1839	3.50	30	244	3.50 × 30
8.007.1840	3.50	35	285	3.50 × 35
8.007.1841	4.00	10	80	4.00 × 10
8.007.1842	4.00	15	121	4.00 × 15
8.007.1843	4.00	20	162	4.00 × 20
8.007.1844	4.00	25	203	4.00 × 25
8.007.1845	4.00	30	244	4.00 × 30
8.007.1846	4.00	35	285	4.00 × 35

2. INTENDED USE / INDICATIONS FOR USE

The HT Supreme™ Drug Coated Coronary Stent System is indicated for improving coronary luminal diameter in patients with symptomatic heart disease due to de novo native coronary artery lesions (length ≤ 35 mm) with reference vessel diameters of 2.25 mm to 4.00 mm.

3. CONTRAINDICATIONS

The HT Supreme™ stent system is contraindicated for use in the following:

- Patients for whom antiplatelet and/or anticoagulant therapy is contraindicated.
- Patients who have lesions that prevent the complete inflation of an angioplasty balloon and/or proper placement of the stent or stent delivery system.
- Patients with known hypersensitivity or contraindication to sirolimus or structurally-related compounds, cobalt, chromium, nickel, tungsten, acrylic-polymers, PLGA polymers or contrast medium.

4. WARNINGS

- The inner packaging must not be opened or damaged prior to use in order to maintain sterility. Use of non-sterile devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.
- Judicious patient selection is necessary because the use of this device has been associated with certain complications, including stent thrombosis, vascular complications, and/or bleeding events. See Section 8 (Potential Adverse Events/Complications) for a full list.
- This product requires the use of anticoagulation and/or antiplatelet therapy. This product should not be used in patients who are not likely to comply with the recommended antiplatelet regimen. Please refer to Section 5.2 (Pre- and Post-Procedure Antiplatelet Regimen) for important information.

5. PRECAUTIONS

5.1 General Precautions

- Stent implantation should only be performed by physicians who have received adequate training.
- Stent implantation should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed in the event of a potential injury or life-threatening complication.
- In cases where multiple stents are implanted, it is recommended that the stent materials be of similar composition. Contact between stents of different materials could lead to increased corrosion potential and/or potential drug interactions in the case of use with other drug eluting stents. Potential drug interactions between the HT Supreme™ stent and other drug eluting stents have not been characterized and should be avoided.
- Subsequent stent blockage or stenosis may require repeat

dilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of endothelialized stent is currently unknown.

- Consideration should be given to the risks and benefits of use in patients with a history of severe reaction to contrast agents.
- Do not expose the delivery system to organic solvents such as alcohol or detergents as the drug carrier coating may be susceptible to damage or premature drug elution.
- Care should be taken to control the guiding catheter tip during stent delivery, deployment, and balloon withdrawal. Before withdrawing the stent delivery system, visually confirm complete balloon deflation by fluoroscopy to avoid guiding catheter movement into the vessel and subsequent arterial damage.
- Stent thrombosis is a rare event but when it occurs it is frequently associated with myocardial infarction (MI) or death. Clinical trials of drug eluting stents evaluated to date have not fully characterized the incidence of stent thrombosis in relation to the risk of cardiac death, MI, or all-cause mortality.
- This device is intended to be used within the specified Indication for Use. Using this device to treat patients and/or lesions outside of the labeled indications may result in an increased risk of adverse events, including stent thrombosis, stent embolization, MI or death.
- Sirolimus is metabolized by the CYP 3A family of enzymes. Studies have shown that drugs such as cyclosporine that inhibit CYP 3A enzymes can slow sirolimus clearance. Use of other "limus" compounds (e.g. everolimus, tacrolimus, etc.) concomitant with a sirolimus eluting stent may compete for binding to metabolic enzymes and have the potential to increase sirolimus levels and reduce clearance rates.
- Use of other drugs that compete for binding to the FKBP12 receptor (e.g. everolimus, tacrolimus, etc.) have the potential to reduce sirolimus efficacy.

5.2 Pre- and Post-Procedure Antiplatelet Regimen

The optimal duration of antiplatelet therapy, specifically P2Y₁₂ inhibitor therapy, is unknown and may not eliminate the risk of thrombosis despite continued therapy. According to recent recommendations provided by the 2016 ACC/AHA Guideline Focused Update on Duration of Dual Antiplatelet Therapy in Patients With Coronary Artery Disease (CAD), continued treatment with a combination of aspirin (81mg) and a P2Y₁₂ inhibitor following stent implantation appears to reduce major adverse cardiac events. Based on clinical evidence including randomized clinical trials, secondary prevention measures, and expert consensus opinion, the guidelines recommend that patients with stable ischemic heart disease undergo aspirin and P2Y₁₂ for at least 6 months in the absence of contraindications; discontinuation after 3 months may be reasonable in patients with a high risk of bleeding. Patients with acute coronary syndrome should receive this treatment for at least 12 months in the absence of contraindications; discontinuation after 6 months may be reasonable in patients with a high risk of bleeding. The complete 2016 ACC/AHA guidelines can be found at: <http://www.onlinejacc.org/content/68/10/1082>

Patient compliance with post-procedural antiplatelet recommendations is very important, as premature discontinuation could result in a higher risk of thrombosis, MI or death. When planning a stent implantation procedure, any future anticipated surgical or dental procedure that could require early discontinuation of antiplatelet therapy should be carefully considered, including whether stent implantation and its associated recommended antiplatelet therapy is the appropriate choice. Conversely, if a surgical or dental procedure that requires suspension of antiplatelet therapy is considered after stent implantation, the risks and benefits of the procedure should be weighed against the possible risks associated with premature discontinuation of antiplatelet therapy. In general, it is recommended to postpone elective surgery for one year or, in cases where surgery cannot be deferred, antiplatelet therapy should be considered during the perioperative period for high risk patients.

In cases where significant active bleeding necessitates the need for early discontinuation of antiplatelet therapy, patients should be monitored carefully for cardiac events and, once stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physicians.

5.3 Use of Multiple Stents

The lesions are to be treated with no more than 1 stent (unless additional stenting is required for bailout situations). In cases where multiple stents are implanted, it is recommended that the stent materials be of similar composition. Contact between stents of different materials could lead to increased corrosion potential due to

the presence of dissimilar metals in a conducting medium. Potential interactions of the HT Supreme™ stent with other drug-eluting or coated stents have not been evaluated and should be avoided whenever possible.

5.4 Intravascular Brachytherapy

The safety and effectiveness of using the HT Supreme™ stent system to treat patients with prior target lesions or patients with prior brachytherapy treatment of previously stented lesions has not been evaluated. Conversely, the effects of brachytherapy treatment on the safety and effectiveness of lesions previously treated with the HT Supreme™ stent system also has not been evaluated. Both treatments may alter arterial remodeling; the interactions between these two treatment has not been established.

5.5 Use In Conjunction With Other Procedures

The safety and effectiveness of using mechanical atherectomy devices (such as directional atherectomy catheters and rotational atherectomy catheters), or laser angioplasty catheters in conjunction with HT Supreme™ stent implantation have not been established.

5.6 Use In Special Populations

5.6.1 Pregnancy

Pregnancy "Category C": See Section 6.6 (Pregnancy) for additional drug information. The HT Supreme™ stent system has not been evaluated in pregnant women. Effects on the developing fetus also have not been evaluated. Effective contraception should be initiated before implanting a HT Supreme™ stent and continued for one year after implantation. While there is no contraindication, the risks and reproductive effects are currently unknown.

5.6.2 Lactation

See Section 6.7 (Lactation) for additional drug information. The effects of the HT Supreme™ stent system on lactation have not been evaluated. A decision should be made whether to discontinue nursing or to implant the stent, taking into account the importance of the stent to the mother.

5.6.3 Ethnicity

Clinical studies of the HT Supreme™ stent system have not specifically evaluated potential differences in safety and effectiveness due to ethnicity.

5.6.4 Pediatric Use

The safety and effectiveness of the HT Supreme™ stent in pediatric patients has not been established.

5.6.5 Geriatric Use

At this time, there is insufficient clinical data to allow differences in the HT Supreme™ stent safety and effectiveness between patients under 65 and over 65 years of age.

5.7 Lesion/Vessel Characteristics

The safety and effectiveness of the HT Supreme™ stent system have not been established for patients with the following target lesion or vessel characteristics or conditions:

- Unresolved vessel thrombus at the lesion site
- Coronary artery reference vessel diameter < 2.25 mm or > 4.00 mm
- Lesion greater than 35mm in length
- Lesions located in saphenous vein grafts, an unprotected left main coronary artery, ostial lesions, chronic total occlusions, and lesions located at a bifurcation
- Previously stented lesions
- Diffuse disease or poor flow (TIMI < 1) distal to the identified lesions
- Excessive tortuosity proximal to or within the lesion
- Recent acute myocardial infarction (AMI) or evidence of thrombus in the target vessel
- Moderate or severe lesion calcification
- Multi-vessel disease
- In-stent stenosis

5.8 Drug Interactions

See Section 6 (Drug Information) for additional information on sirolimus drug interactions. The effects of potential drug interactions on the safety and effectiveness of the HT Supreme™ stent system have not been evaluated. The effectiveness of sirolimus is based on its binding to the FKBP12 receptor and subsequent interaction with the mechanistic Target of Rapamycin (mTOR). Concomitant use of other drugs that bind this receptor can potentially reduce sirolimus efficacy. Several drugs are known to affect the metabolism of sirolimus, and other drug interactions may be inferred from known

metabolic effects. Sirolimus is known to be a substrate for both cytochrome P450 IIIA4 (CYP3A4) and P-glycoprotein. Consideration should be given to the potential effect of drug interaction when deciding to implant a HT Supreme™ stent in a patient who is taking a drug that could interact with sirolimus by competing for sirolimus metabolism or receptor binding, or when deciding to initiate a therapy with such a drug in a patient who had been implanted with HT Supreme™ stent.

5.9 Immune Suppression Potential

Sirolimus, the HT Supreme™ stent active ingredient, is an immunosuppressive agent. Potential suppression of the immune system in patients treated with the HT Supreme™ stent system has not been specifically evaluated. Systemic immune suppression is unlikely with single stent implantation due to the small quantity of sirolimus on a single stent. If, however, patients are treated with multiple HT Supreme™ stent systems, it may be possible for sirolimus systemic concentrations to temporarily approach immunosuppressive levels, especially in patients who also have hepatic insufficiency or who are taking drugs that inhibit CYP3A4 or P-glycoprotein. For this reason, consideration should be given to patients taking other immunosuppressive agents or who are at risk for immune suppression.

5.10 Lipid Elevation Potential

The use of oral sirolimus in renal transplant patients was associated with increased serum cholesterol and triglycerides that in some cases required treatment. The effect was seen with both low and high dose prolonged oral therapy in a dose related manner. When used according to the indications for use, the systemic sirolimus concentrations from the HT Supreme™ stent are expected to be lower than the concentrations usually obtained in transplant patients, but the magnitude and duration of any effect of those concentrations on lipids is not known.

5.11 Magnetic Resonance Imaging (MRI)

In non-clinical testing, the HT Supreme™ stent has been shown to be MR Conditional for single stents up to 35 mm in length or overlapped stent lengths up to 67 mm. A patient with this device can be safely scanned in an MR system that meets the following conditions:

- Static magnetic field of 1.5 T and 3.0 T;
- Maximum spatial field gradient of 4,000 gauss/cm or less;
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of ≤ 2 W/kg (Normal Operating Mode) for a 1.5 T or a 3.0 T MRI system.

Under the scan conditions defined above, non-clinical testing results (including the effect of blood flow) indicate the HT Supreme™ stent is expected to produce a maximum temperature rise of less than 5.0°C after 15 minutes of continuous scanning. In non-clinical testing, the image artifact caused by the device extends approximately 6 mm from the HT Supreme™ stent system when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

5.12 Stent Handling

- The HT Supreme™ stent is for single use only. Do not resterilize or reuse this device. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. SINOMED will not be responsible for any direct, incidental or consequential damages resulting from resterilization reuse.
- Note the "Use by" (expiration) date on the product label.
- The outside surface of the foil pouch is NOT sterile. Only the contents inside the foil pouch should be considered sterile.
- Do not remove the stent from the delivery system. Removal may damage the stent, the stent coating, and/or lead to stent embolization. These components are intended to perform together as a single system. Do not use the stent with a different delivery system.
- Do not use the delivery system in conjunction with other stents.
- Special care must be taken not to handle or disrupt the stent on the balloon especially during delivery system removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub as it may damage the stent, polymer system and/or lead to stent embolization.
- Do not manipulate, touch, or handle the stent with your fingers, which may cause coating damage, contamination, or stent dislodgement from the delivery balloon.

- Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in stent deployment. See Section 12.3.3 (Operator's Instructions; Delivery System Preparation).

5.13 Stent Placement

5.13.1 Stent Delivery System Preparation

- Do not prepare or pre-inflate the delivery system prior to stent deployment other than as directed. This may cause dislodgement of the stent from the balloon. Use the balloon purging technique described in Section 12.3.3 (Operator's Instructions; Delivery System Preparation).
- While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.
- Use guiding catheters which have lumen sizes that are suitable to accommodate the stent delivery system. See Section 1.1 (Device Component Description).

5.13.2 Stent Implantation

- Pre-dilate the target vessel with an appropriately sized angioplasty balloon. Failure to do so may increase the difficulty of stent placement and cause procedural complications.
- Do not expand the stent if it is not properly positioned in the vessel. See Section 5.13.3 (Stent System Removal)
- Implanting a stent may lead to vessel dissection and acute closure requiring additional intervention (CABG, further dilatation, placement or additional stents, or other).
- The safety and effectiveness of treating more than one vessel per coronary artery with HT Supreme™ stents has not been established. In cases where this is performed, it is recommended that the stent be placed in the distal lesion first in order to minimize dislodgement risk incurred by traversing through deployed stents.
- Stent placement may compromise side branch patency.
- Do not exceed Rated Burst Pressure (RBP) as indicated on product label. See Section 13 (In Vitro Compliance Information) for stent compliance information as a function of diameter. Use of a pressure monitoring device is recommended to prevent over pressurization. Applying pressures higher than specified on the product label may result in a ruptured balloon with possible arterial damage and dissection. The stent inner diameter should approximate 1.1 times the reference diameter of the vessel. The size of stent expansion should be confirmed by radiography.
- An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed as the stent may be damaged when retracting the undeployed stent back into the guiding catheter.
- Should any resistance be felt at any time during stent system withdrawal, the stent delivery system and guiding catheter should be removed as a single unit. See Section 5.13.3 (Stent System Removal).
- Although the stent delivery system balloon is strong enough to expand the stent without rupture, a circumferential balloon tear distal to the stent and prior to complete stent expansion, could cause the balloon to become tethered to the stent, requiring surgical removal. In case of balloon rupture, it should be withdrawn and, if necessary, a new dilatation catheter exchanged over the guide wire to complete the expansion of the stent.
- Ensure the stented area covers the entire lesion/ dissection site and that no gaps exist between stents.

5.13.3 Stent System Removal

- Should any resistance be felt at any time during either lesion access or when removing the delivery system post-stent implantation, the stent delivery system and the guiding catheter should be removed as a single unit by performing the following steps under direct fluoroscopic visualization:
 - Confirm complete balloon deflation. The deflation time is no more than 20 seconds for all sizes. If unusual resistance is felt during delivery system withdrawal, pay particular attention to the guiding catheter position. In some cases, it may be necessary to slightly retract the guiding catheter in order to prevent unplanned guiding catheter movement, a coronary tree angiographic assessment should be undertaken to ensure that there is no damage to the coronary vasculature.
 - DO NOT retract the delivery system into the guiding catheter.

- Position the proximal balloon marker just distal to the guiding catheter tip.
- Advance the guide wire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating hemostatic valve to secure the delivery system to the guiding catheter, and remove the guiding catheter and delivery system as a single unit.
- Failure to follow these steps and/or applying excessive force to the delivery system can potentially result in loss or damage to the stent and/or delivery system components and/or vasculature. If it is necessary to retain guide wire position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components. Stent retrieval methods (i.e., additional wires, snares and/or forceps) may result in additional trauma to the artery vasculature and/or the vascular access site. Complications may include, but are not limit to, bleeding, hematoma, or pseudoaneurysm.

5.13.4 Post-Procedure

- When crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guide wire, a balloon catheter or delivery system, exercise care to avoid disrupting the stent placement, apposition, geometry and/or coating.
- Antiplatelet therapy should be administered post-procedure. See Section 5.2 (Pre- and Post-Procedure Antiplatelet Regimen). At risk patients who require early discontinuation of antiplatelet therapy should be monitored carefully for cardiac events. At the discretion of the patient's treating physician, the antiplatelet therapy should be restarted as soon as possible.
- Post-procedure, patients can safely undergo an MRI under certain conditions. Please refer to Section 5.11 (Magnetic Resonance Imaging (MRI)) if a post-procedure MRI is needed.

6. DRUG INFORMATION

6.1 Mechanism of Action

Sirolimus is an anti-proliferative drug that binds to FK506-binding protein 12 (FKBP12). The rapamycin-FKBP12 dimeric complex inhibits the mechanistic target of rapamycin (mTOR). mTOR is one of the key regulatory proteins that directs cell metabolism and cell growth, in part, by upregulating the cyclin-dependent kinase inhibitor p27. Upregulation of p27 leads to inhibition of cell cycle progression at the juncture of G1 to S phase. Thus, in most cells, therapeutic concentrations of sirolimus will arrest cell division at the late G1 stage. Because this is prior to initiation of mitosis, the arrest at G1 causes a cytostatic rather than cytotoxic inhibition of cell growth. The anti-proliferative properties of sirolimus make it a useful agent to reduce the risk of restenosis. In this capacity, it interferes with fibroblast and smooth muscle cell proliferation as well as the proliferation of inflammatory cells. These three cell types are largely responsible for the neointimal hyperplasia resulting from vascular injury secondary to stent implantation.

6.2 Pharmacokinetics

Sirolimus Pharmacokinetics (PK) when eluted from the HT Supreme™ stent post-implantation has been evaluated in animals. Whole blood sirolimus PK results determined in the porcine animal models are provided in Table 6. 1. BQL referred to in the table stands for 'below the limit of quantification.'

Table 6.1 Whole Blood Pharmacokinetics Results for HT Supreme™ (ng/ml)

0 min	5 min	15 min	30 min	60 min	180 min
BLQ	1.35 ± 0.27	1.49 ± 0.23	1.39 ± 0.22	1.36 ± 0.23	1.39 ± 0.20

1 Day	3 Days	7 Days	14 Days
0.156 ± 0.087	0.038 ± 0.059	<0.1	0.036 ± 0.056

28 Days	60 Days	90 Days	120 Days	150 Days
0.094 ± 0.076	BLQ	BLQ	BLQ	BLQ

Note: Blood drug levels were determined from an "n" of 6; BLQ < 0.1 ng/ml

The drug levels found in whole blood after stent implantation demonstrated a rapid increase at 5 minutes, with the concentration maintained at the same level for at least 180 minutes (1.35-1.39 ng/mL). This was followed by a decrease to 0.156 ng/mL at 1 day and to below detection level (0.1 ng/mL) at 60 days. Overall, there was a 90% decrease in sirolimus blood concentration at 1 day in comparison to the first three hours post stent implantation. During this PK study, the maximum observed blood concentration of sirolimus never reached the minimum therapeutic value of 3.0 ng/mL necessary for effective systemic administration to prevent organ rejection in patients. Sirolimus disappearance from

circulation following stent implantation should further limit systemic exposure and adverse events associated with long-term systemic administration at therapeutic levels. Despite limited systemic exposure to sirolimus, consistent local arterial delivery of sirolimus from the stent has been demonstrated in pre-clinical studies.

6.3 Interactions with Drugs or other Substances

Formal drug interaction studies have not been performed with the HT Supreme™ stent because of limited systemic exposure to sirolimus following implantation of the stent. See Section 5.8 (Drug Interactions) and Section 6.2 (Pharmacokinetics). However, consideration should be given to the potential for both systemic and local drug interactions in the vessel wall when deciding to place the HT Supreme™ stent in a patient taking a drug with known interaction with sirolimus. The amount of drug that circulates in the bloodstream following implantation of a HT Supreme™ stent is significantly lower than that obtained with oral doses. Sirolimus is extensively metabolized by cytochrome P450 3A4 (CYP3A4) in the gut wall and liver and undergoes efflux from enterocytes of the small intestine by P-glycoprotein (P-gp). Therefore, absorption and the subsequent elimination of systemically absorbed sirolimus may be influenced by drugs that affect these proteins. Inhibitors of CYP3A4 and P-gp may increase sirolimus levels, while inducers of CYP3A4 and P-gp may decrease sirolimus levels. The pharmacokinetic interaction between orally administered sirolimus and concomitantly administered drugs is discussed below. Drug interaction studies have not been conducted with drugs other than those described below.

6.4 Ketoconazole

Multiple-dose ketoconazole administration significantly affected the rate and extent of absorption and sirolimus exposure after administration of a sirolimus oral formulation, as reflected by increases in sirolimus C_{max} , t_{max} , and AUC of 4.3-fold, 38%, and 10.9-fold, respectively. However, the terminal $t_{1/2}$ of sirolimus was not changed. Single-dose sirolimus did not affect steady-state 12-hour plasma ketoconazole concentrations. It is recommended that sirolimus oral solution and oral tablets should not be administered with ketoconazole.

6.4.1 Rifampin

Pre-treatment of 14 healthy volunteers with multiple doses of rifampin, 600 mg daily for 14 days, followed by a single 20-mg dose of sirolimus, greatly increased sirolimus oral-dose clearance by 5.5-fold (range = 2.8 to 10), which represents mean decreases in AUC and C_{max} of about 82% and 71%, respectively. In patients where rifampin is indicated, alternative therapeutic agents with less enzyme induction potential should be considered.

6.4.2 Diltiazem

The simultaneous oral administration of 10 mg of a sirolimus oral solution and 120 mg of diltiazem to 18 healthy volunteers significantly affected the bioavailability of sirolimus. Sirolimus C_{max} , t_{max} , and AUC were increased 1.4-, 1.3-, and 1.6-fold, respectively. Sirolimus did not affect the pharmacokinetics of either diltiazem or its metabolites desacetyldiltiazem and desmethylidiltiazem.

6.4.3 Cyclosporine

Single-dose pharmacokinetic interactions between cyclosporine and sirolimus were investigated for two sirolimus oral formulations in studies using 24 healthy volunteers. Compared to results obtained when oral sirolimus was administered alone, the oral administration of 10 mg sirolimus 4 hours after a single dose of 300 mg cyclosporine soft gelatin capsules increased mean sirolimus AUC by 33% to 80% and increased mean sirolimus C_{max} by 33% to 58%, depending on the sirolimus formulation. The half-life of sirolimus was not significantly affected. The cyclosporine mean AUC and mean C_{max} were not significantly affected.

6.4.4 Drugs which may be coadministered without dose adjustment

Clinically significant pharmacokinetic drug-drug interactions were not observed in studies of drugs listed below in conjunction with orally administered sirolimus. Sirolimus and these drugs may be coadministered without dose adjustments.

- Acyclovir
- Digoxin
- Glyburide
- Nifedipine
- Norgestrel/ethinyl estradiol
- Prednisolone
- Sulfamethoxazole/trimethoprim

6.4.5 Other drug interactions

Drugs that may increase sirolimus blood concentrations include:

- Calcium channel blockers: nicardipine, verapamil.
- Antifungal agents: clotrimazole, fluconazole, itraconazole.
- Macrolide antibiotics: clarithromycin, erythromycin, troleandomycin.
- Gastrointestinal prokinetic agents: cisapride, metoclopramide.
- Other drugs: bromocriptine, cimetidine, danazol, HIV-protease inhibitors (e.g., ritonavir, indinavir).
- Drugs that may decrease sirolimus levels include:
- Anticonvulsants: carbamazepine, phenobarbital, phenytoin.
- Antibiotics: rifabutin, rifampentine.

These lists are not all inclusive.

Care should be exercised when drugs or other substances that are metabolized by CYP3A4 are administered concomitantly with implantation of the HT Supreme™ stent.

6.4.6 Grapefruit juice

Grapefruit juice reduces CYP3A4-mediated metabolism of sirolimus.

6.4.7 Herbal Preparations

St. John's Wort (*Hypericum perforatum*) induces CYP3A4 and P-glycoprotein. Because sirolimus is a substrate for both cytochrome CYP3A4 and P-glycoprotein, there is the potential that the use of St. John's Wort in patients receiving the HT Supreme™ stent could result in reduced sirolimus levels.

6.4.8 Vaccination

Immunosuppressants may affect response to vaccination. Therefore, for some period after receiving a HT Supreme™ stent, vaccination may be less effective. The use of live vaccines should be avoided; live vaccines may include, but are not limited to, measles, mumps, rubella, oral polio, BCG, yellow fever, varicella, and TY21a typhoid.

6.4.9 Drug-laboratory test interactions

There are no studies on the interactions of sirolimus in commonly employed clinical laboratory tests.

6.5 Carcinogenicity, Genotoxicity, and Reproductive Toxicity

All of the components of the HT Supreme™ stent have been used previously on other coronary drug-eluting stents. The cobalt chromium alloy used as the stent platform is commonly used in most coronary artery stents. The HT Supreme™ stent has a similar amount of sirolimus as the previously commercialized Cypher® sirolimus eluting stent. The base coat poly butyl methacrylate (PBMA) polymer is the same as that used on the Cypher® and also the Xience®, Promus™, and Endeavor Resolute® drug eluting stents. The top coat biodegradable poly lactide-co-glycolic acid (PLGA) polymer is also used on the Synergy® stent.

Therefore, previous testing that was performed on these stents is applicable to the HT Supreme™ stent.

The HT Supreme™ stent has been evaluated for genotoxicity in three different assay systems: The Ames bacterial reverse mutation assay, the in vitro mouse lymphoma assay and the in vivo mouse micronucleus assay. Results of these assays show the HT Supreme™ stent is not mutagenic nor clastogenic in nature.

The sirolimus component of the HT Supreme™ stent has been evaluated for carcinogenic and reproductive toxicity potential with respect to its use in much higher doses as an immunosuppressant. These data are available as part of the Rapamune® product testing. Carcinogenicity studies were conducted in female mice and male and female rats. In the 86-week female mouse study at dosages of 0, 12.5, 25 and 50/6 mg/kg/day (dosage lowered from 50 to 6 mg/kg/day at week 31 due to infection secondary to immunosuppression) there was a statistically significant increase in malignant lymphoma at all dosages (approximately 86 to 357 times the maximum recommended human dose [MRHD]) compared to controls. In the 104-week rat study at dosages of 0, 0.05, 0.1, and 0.2 mg/kg/day, there was an increased incidence of testicular adenoma in the 0.1 and 0.2 mg/kg/day (approximately 1.4 to 2.9 times the MRHD) groups.

There was no effect on fertility in female rats following the administration of sirolimus at dosages up to 0.5 mg/kg (approximately 7 times the MRHD). In male rats, there was no significant difference in fertility rate compared to controls at a dosage of 2 mg/kg (approximately 28 times the MRHD). Reductions in testicular weights and/or histological lesions (e.g., tubular atrophy and tubular giant cells) were observed in rats following dosages of 0.65 mg/kg (approximately 9 times the MRHD) and above and in a monkey study at 0.1 mg/kg (approximately 1.4 times the MRHD)

and above. Sperm counts were reduced in male rats following the administration of sirolimus for 13 weeks at a dosage of 6 mg/kg (approximately 85 times the MRHD), but showed improvement by 3 months after dosing was stopped.

Based on these data, and the fact that much higher doses of sirolimus have been safely administered systemically as an immunosuppressant, the total drug load of the HT Supreme™ stent would not be expected to result in adverse long term effects including carcinogenicity. Neither the base coat nor the top coat polymer have shown any indication of genotoxicity, carcinogenicity or reproductive toxicity based on published literature. Their previous history of safe clinical use as coatings for other drug eluting stents suggests that they would not contribute to long term genotoxic effects. The PBMA polymer is a permanent coating and the PLGA polymer degrades to lactic acid and glycolic acid, two naturally occurring compounds that are readily metabolized to carbon dioxide and water by normal metabolic pathways.

6.6 Pregnancy

Pregnancy "Category C": There are no adequate and well controlled studies in pregnant women of sirolimus or HT Supreme™ stents. In rabbit studies, there was no effect of sirolimus on development at the maternally toxic dosage of 0.05 mg/kg/day (approximately 3 times the dose provided by a stent coated with 285µg sirolimus adjusted for body surface area). Effective contraception should be initiated before implanting a HT Supreme™ stent and continued for one year after implantation. The HT Supreme™ stent should be used during pregnancy only if the potential benefit outweighs the potential risks.

6.7 Lactation

It is unknown whether sirolimus is distributed in human milk. Also, sirolimus pharmacokinetic and safety profiles have not been determined in infants. Consequently, mothers should be advised of potential serious adverse reactions to sirolimus in nursing infants. Prior to HT Supreme™ stent implantation, decisions should be made regarding whether to discontinue nursing or conduct an alternate percutaneous coronary intervention procedure.

7. OVERVIEW OF CLINICAL STUDIES

Currently there are one completed trial and three on-going randomized prospective trials using the 2nd generation HT Supreme™ stent. The PIONEER I study in Europe is fully enrolled (170 subjects) and 3 years follow up data was available. The PIONEER II, PIONEER II-OCT and PIONEER III studies are being conducted worldwide.

7.1 The PIONEER I Study

The PIONEER I study is a European-based prospective, multi-center, randomized First-in Man trial to compare the safety and efficacy of the HT Supreme™ DES system biodegradable polymer sirolimus-eluting stent against the Resolute Integrity zotarolimus-durable polymer stent in patients with de novo coronary artery lesions. The PIONEER I clinical trial, randomized 1:1, in a single blind fashion, is being conducted at 14 interventional cardiology centers in The Netherlands, Belgium, Spain and Portugal. Clinical follow-up occurred at 1, 9, 24 and 36 months' post-stent implantation. All subjects' angiography (QCA) assessment was performed at baseline (pre- and post-procedure) and at the 9 months' follow-up visit. Subjects ("n" = 170); having documented stable or unstable angina or silent ischemia) were enrolled with a maximum of two de novo native coronary artery lesions located in two different epicardial vessels. The study's primary endpoint is in-stent Late Lumen Loss (LLL) at 9 months after stent implantation as assessed by off-line QCA.

The study is now considered complete with regard to the primary endpoints and 3 years follow-up.

7.2 The PIONEER II Study

The PIONEER II study is a Prospective multi-center clinical trial being conducted in China to evaluate the efficacy and safety of the HT Supreme™ DES biodegradable polymer sirolimus-eluting stent in patients with de novo coronary lesions. The study is comprised of 2 arms: The first arm is a randomized controlled trial comparing the HT Supreme™ DES to the first generation BUMA as the control article.

The primary end-point is 9-month QCA driven late lumen loss. The second arm is designed to generate Objective Performance Criteria (OPC) with a primary endpoint of target lesion failure at 1 year. All subjects enrolled will have clinical follow-up out to 5 years. Currently the PIONEER II study completed the patients enrollment in both the RCT and OPC arms in 2018, and the 9-month QCA angiographic follow up and 1 year follow up is in process. There is no data available at this time.

7.3 The PIONEER II OCT Study

PIONEER II OCT is a prospective, multicenter, randomized study, comparing the HT Supreme stent versus the Xience stent on the extent of neointima formation at 1- and 2- months using Optical Coherence Tomography (OCT) imaging after stent implantation in high bleeding risk CAD patients.

The study aimed to provide the OCT imaging data of the HT Supreme stent compared to the durable polymer coated Xience family of stents, to provide evidence of early healing capacity of the HT Supreme stent based on OCT imaging technique. The primary endpoint was the rate of strut coverage at 1- or 2- month.

A total of 82 patients were enrolled in the PIONEER II OCT study at 8 sites in China. A total of 75 patients were available for OCT follow-up: 16 patients (including 18 lesions) were randomized to treatment with the HT Supreme stent and 15 patients (including 17 lesions) were randomized to treatment with the Xience stent in cohort-1 with a 1 month OCT analysis; 21 patients (including 21 lesions) were randomized to treatment with the HT Supreme stent and 23 patients (including 28 lesions) were randomized to treatment with the Xience stent in cohort-2 with a 2 month OCT analysis.

One-month (Cohort-1) OCT follow-up results demonstrates that the percentage of strut coverage was significantly higher in the HT Supreme than the Xience group (83.8% vs. 73.0%, $P_{\text{superiority}} = 0.037$), while the mean neointimal area was greater in the HT Supreme group compared to the Xience group (0.63 mm² vs. 0.36 mm², $P=0.008$).

Two-months (Cohort-2) OCT follow-up results demonstrates that the percentage of strut coverage was higher with the HT Supreme stent than the Xience group (80.3% vs. 73.3%, $P_{\text{noninferiority}}=0.006$), while the mean neointimal area was greater in the HT Supreme group compared to the Xience group (0.63 mm² vs. 0.49 mm², $P=0.16$).

7.4 The PIONEER III Study

The PIONEER III is a prospective multi-center clinical trial being conducted in the United States, Canada, Japan and Europe. The PIONEER III study is to evaluate the efficacy and safety of the HT Supreme stent in patients with de novo coronary lesions comparing with the Xience/Promus family stents. The primary analysis will be a non-inferiority test comparing HT Supreme DES to Durable Polymer coated Everolimus Eluting Stent (DP EES) for the primary safety and efficacy endpoint of target lesion failure at 12 months. A total of 1632 patients (2:1 randomization of HT Supreme vs. the control group) will be enrolled at up to 100 sites worldwide.

Currently the PIONEER III study completed the patients enrollment in 2019, and 1 year follow up is in process. There is no data available as of now.

8. POTENTIAL ADVERSE EVENTS/ COMPLICATIONS

The potential adverse events or complications (listed in alphabetical order) which may be associated with the use of coronary stenting devices, IVUS (intravascular ultrasound) or PCI (percutaneous coronary intervention) include, but are not limited to:

- Allergic reaction (to contrast, antiplatelet therapy, stent system-material, drug, or polymer coating);
- Arrhythmias;
- Aneurysm;
- Coronary artery spasm;
- Death;
- Dissection;
- Failure to deliver the stent;
- Hemorrhage requiring transfusion;
- Hypotension/Hypertension;
- Incomplete stent apposition;
- Infection and pain at the vascular access site;
- Stent embolization.

The occurrence of the above listed complications may lead to the need for a repeat catheterization and/or percutaneous coronary intervention, myocardial infarction, emergency bypass surgery, or death. The following additional side effects/complications (in alphabetical order) may be associated with, but not limited to, the use of sirolimus:

- Abnormal liver function tests;
- Anemia;
- Diarrhea;
- Hypercholesterolemia;
- Hypersensitivity, including anaphylactic / anaphylactoid type reactions;

- Hypertriglyceridemia;
- Infections;
- Leucopenia;
- Renal failure;
- Pseudoaneurysm;
- Thrombocytopenia.

The side effects/complications of the HT Supreme™ stent polymer are no different than those of other stent coatings and may include, but are not limited to, the following (in alphabetical order):

- Allergic reaction
- Focal inflammation at the site of stent implantation
- Restenosis of the stented artery
- There may be other adverse events or complications that are unforeseen at this time.

9. INDIVIDUALIZATION OF TREATMENT

See also Section 5 (Precautions). The risks and benefits should be carefully considered for each patient before use of the HT Supreme™ stent system. Among patient selection factors to be assessed includes a judgment regarding risk of prolonged anticoagulation. Stenting should be generally avoided in those patients at heightened risk of bleeding (e.g., those patients with recently active gastritis or peptic ulcer disease). Antiplatelet drugs should be used in combination with the HT Supreme™ stent. Physicians should adhere to current ACC/AHA/SCAI Guidelines as described in Section 5.2 (Pre- and Post-Procedure Antiplatelet Regimen). The specific needs of individual patients should be considered to determine the specific antiplatelet/anticoagulation regimen to be used for patients in general practice.

Each HT Supreme™ stent must be sized appropriately to fit the patient's particular coronary anatomy. Proper sizing of the device (diameter and length) is the responsibility of the physician.

The risks and benefits previously described should be carefully considered for each patient before use of the HT Supreme™ stent system. Co-morbidities that increase the risk of poor initial results in PCI or increase the risks of emergency referral for bypass surgery should be reviewed.

10. PATIENT COUNSELING AND PATIENT INFORMATION

Physicians should consider the following in counseling patients about this product:

- Discuss the risks associated with stent placement.
- Discuss the risks associated with a sirolimus-eluting stent.
- Discuss the risks of early discontinuation of the antiplatelet therapy.
- Discuss the risks of late stent thrombosis with drug eluting stent use in higher risk patient subgroups.
- Discuss the risk/benefit issues for this particular patient.
- Discuss alteration to current life-style immediately following the procedure and over the long term.

The following patient materials are available for this product:

- A Patient Information Guide, which includes information on coronary artery disease, the implant procedure and the HT Supreme™ stent system (provided to physician)
- A Stent Implant Card that includes both patient information and stent implant information (provided in package)

11. HOW SUPPLIED

- Sterile: This device is sterilized with E-Beam irradiation. Non-pyrogenic. It is intended for single use only. Do not resterilize. Do not use if the package is opened or damaged.
- Contents: One (1) HT Supreme™ stent system and one (1) Stent Implant Card.
- Storage: Store in the original container in a dry, dark, cool place. Store between 10°C and 30°C.
- Disposal: After use, this product and its packaging should be disposed of in accordance with standard policies and procedures required by the hospital, administration, and/or local government.

12. OPERATOR'S INSTRUCTIONS

Prepare the guiding catheter and guidewire according to the manufacturer's instructions. The HT Supreme™ stent system is compatible with 0.36 mm (0.014 in) guidewires. Refer to product labeling or Section 12.2 (Materials Required) for specific guiding catheter compatibility.

12.1 Inspection Prior to Use

- Carefully inspect the sterile package before opening and check for damage to the sterile barrier. Do not use if the integrity of the sterile package has been compromised.
- Tear open the foil pouch and remove the inner pouch. Note: The outside surface of the foil pouch is NOT sterile. Only the contents of the foil pouch should be considered sterile.
- Prior to using the HT Supreme™ stent system, carefully remove the system from the package and inspect for bends, kinks, and other damage. Verify that the stent does not extend beyond the radiopaque balloon marker. Do not use if any defects are noted. However, do not manipulate, touch, or handle the stent with your fingers, which may cause coating damage, contamination, or stent dislodgement from the delivery balloon.

Note: At any time during use of the HT Supreme™ stent system, if the stainless steel proximal shaft has been bent or kinked, do not continue to use the catheter.

12.2 Material Required

- 1.42 mm (5 Fr/0.056 in) minimum inner diameter guide catheter
- 10 - 20 cc syringe
- Heparinized normal saline
- 0.36 mm (0.014 in) maximum outer diameter guidewire
- Rotating hemostatic valve with appropriate minimum inner diameter (0.096 inch [2.44 mm])
- Contrast medium diluted 1:1 with heparinized normal saline (HepNS)
- Inflation device
- Torque device
- Three-way stopcock
- Pre-deployment dilatation catheter
- Guidewire introducer
- Appropriate arterial sheath
- Appropriate anticoagulation and antiplatelet drugs

12.3 Preparation

12.3.1 Packaging Removal

Note: The outside surface of the foil pouch is NOT sterile. Only the contents of the foil pouch should be considered sterile.

1. Carefully remove the delivery system from its protective tubing for preparation of the delivery system. Do not bend or kink the hypotube during removal.
2. Remove the product mandrel and protective stent sheath by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector and gently remove distally. If unusual resistance is felt during product mandrel and stent sheath removal, do not use this product and replace with another. Follow product returns procedure for the unused device.

12.3.2 Guide Wire Lumen Flush

Flush the guide wire lumen with HepNS until fluid exits the guide wire exit notch.

Note: Avoid manipulation of the stent while flushing the guide wire lumen, as this may disrupt the placement of the stent on the balloon.

12.3.3 Delivery System Preparation

1. Prepare an inflation device/syringe with diluted contrast medium.
2. Attach an inflation device/syringe to the stopcock; attach it to the inflation port of the product. Do not bend the product hypotube when connecting to the inflation device/syringe.
3. With the tip down, orient the delivery system vertically.
4. Open the stopcock to the delivery system, pull negative for 30 seconds, release to neutral for contrast fill.
5. Close the stopcock to the delivery system; purge the inflation device/syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled. If bubbles persist, do not use the product.
7. If a syringe was used, attach a prepared inflation device to stopcock.
8. Open the stopcock to the delivery system.
9. Leave on neutral

Note: While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.

Note: If air is seen in the shaft, repeat Section 12.3.3 (Delivery System Preparation), steps 3 through 5 to prevent uneven stent expansion.

12.3.4 Delivery Procedure

1. Prepare the vascular access site according to standard practice.
 2. Pre-dilate the lesion with a PTCA catheter of appropriate length and diameter for the vessel/lesion to be treated. Limit the longitudinal length of pre-dilatation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the HT Supreme™ stent.
- Note:** The labeled stent diameter refers to expanded stent inner diameter.
3. Maintain neutral pressure on the inflation device attached to the delivery system. Open the rotating hemostatic valve as wide as possible.
 4. Backload the delivery system onto the proximal portion of the guide wire while maintaining guide wire position across the target lesion.
 5. Carefully advance the delivery system into the guiding catheter and over the guide wire to the target lesion. When using a Rapid Exchange (RX) system be sure to keep the hypotube straight. Ensure guiding catheter stability before advancing the stent system into the coronary artery.

Note: If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance may indicate a problem and the use of excessive force may result in stent damage or dislodgement. Maintain guide wire placement across the lesion and remove the delivery system and guiding catheter as a single unit.

6. Advance the delivery system over the guide wire to the target lesion under direct fluoroscopic visualization. Utilize the radiopaque balloon markers to position the stent across the lesion. Perform angiography to confirm stent position. If the position of the stent is not optimal, it should be carefully repositioned or removed (see Section 5.13.3 (Stent System Removal)). The balloon markers indicate both the stent edges and the balloon shoulders. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion.

Note: Should any resistance be felt at any time during either lesion access or removal of the delivery system post stent implantation, remove the entire system as a single unit. See Section 5.13.3 (Stent System Removal) for specific delivery system removal instructions.

7. Tighten the rotating hemostatic valve. The stent is now ready to be deployed.

12.4 Deployment Procedure

CAUTION: Refer to Table 13.1 HT Supreme™ Stent Compliance Inflation Pressure vs. Stent Inner Diameter.

1. Prior to deployment, reconfirm the correct position of the stent relative to the target lesion using the radiopaque balloon markers.
2. Deploy the stent slowly by pressurizing the delivery system in 2 atm increments, every 5 seconds, until stent is completely expanded. Accepted practice generally targets an initial deployment pressure that would achieve a stent inner diameter ratio of about 1.1 times the reference vessel diameter. (See Table 13.1). Maintain pressure for 30 seconds. If necessary, the delivery system can be repressurized or further pressurized to assure complete apposition of the stent to the artery wall. Do not exceed the labeled rated burst pressure (RBP) of 18 atm for 2.25 mm to 3.50 mm diameters and 16 atm for 4.0 mm diameter.
3. Fully cover the entire lesion and balloon treated area (including dissections) with the HT Supreme™ Stent, allowing for adequate stent coverage into healthy tissue proximal and distal to the lesion.
4. Deflate the balloon by pulling negative on the inflation device for 30 seconds. Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during stent delivery system withdrawal, pay particular attention to guiding catheter position.
5. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography or intravascular ultrasound (IVUS).
6. If the deployed stent size is still inadequate with respect to reference vessel diameter, a larger balloon may be used to further expand the stent. If the initial angiographic appearance is sub-optimal, the stent may be further expanded using a low profile, high pressure, non-compliance balloon dilatation catheter. If this is

required, the stented segment should be carefully recrossed with a prolapsed guide wire to avoid disrupting the stent geometry. Deployed stents should not be left underdilated.

- If more than one HT Supreme™ stent is needed to cover the lesion and balloon treated area, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped. To ensure that there are no gaps between stents the balloon marker bands of the second HT Supreme™ stent should be positioned inside the deployed stent prior to expansion.
- Reconfirm stent position and angiographic results. Repeat inflations until optimal stent deployment is achieved.

12.5 Removal Procedure

- Deflate the balloon by pulling negative pressure on the inflation device for 30 seconds. Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during stent delivery system withdrawal, pay particular attention to the guiding catheter position.
- Fully open the rotating hemostatic valve.
- While maintaining the guide wire position and negative pressure on the inflation device, withdraw the delivery system.

Note: Should any resistance be felt at any time during either lesion access or removal of the delivery system post stent implantation, the entire system should be removed as a single unit. See Section 5.13.3 (Stent System Removal) for specific delivery system removal instructions.

- Tighten the rotating hemostatic valve.
- Repeat angiography to assess the stented area. If post-dilatation is necessary, ensure that the final stent diameter matches the reference vessel diameter. Assure that the stent is not underdilated.

12.6 Post-Deployment Dilatation of Stent Segments

All efforts should be taken to assure that the stent is not underdilated. If the deployed stent size is still inadequate with respect to the vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. The stent may be further expanded using a low profile, high pressure, and non-compliant balloon catheter. If this is required, the stented segment should be recrossed carefully with a prolapsed guide wire to avoid dislodging the stent. The balloon should be centered within the stent and should not extend outside of the stented region.

CAUTION: Do not dilate the stent beyond the following maximum stent inner diameter.

Nominal Stent Diameter (mm)	Maximum Stent Inner Diameter (mm)
2.25 to 3.00	3.34 mm
3.25 to 4.00	4.14 mm

13. IN VITRO COMPLIANCE INFORMATION

Table 13.1: HT Supreme™ Stent Compliance Inflation Pressure vs. Stent Inner Diameter

Inflation Pressure		Stent Inner Diameter (mm)							
atm	kPa	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00
8	811	2.14	2.38	2.56	2.81	3.07	3.29	3.59	
9	912	2.20	2.45	2.65	2.90	3.16	3.40	3.71	
10*	1013*	2.25	2.50	2.75	3.00	3.25	3.50	3.80	
11	1115	2.27	2.57	2.78	3.03	3.29	3.56	3.90	
12*	1216*	2.29	2.62	2.83	3.09	3.32	3.63	4.00	
13	1317	2.31	2.66	2.88	3.14	3.37	3.69	4.04	
14	1419	2.33	2.70	2.92	3.18	3.41	3.74	4.06	
15	1520	2.35	2.73	2.96	3.23	3.46	3.79	4.10	
16**	1621**	2.38	2.77	3.00	3.26	3.49	3.85	4.14	
17	1723	2.40	2.81	3.03	3.30	3.53	3.91		
18**	1824**	2.43	2.85	3.07	3.34	3.57	3.98		

Note: This information is based on in vitro testing at 37°C and do not take into account lesion resistance. At least 99.9 percent of balloons (with a 95 percent confidence) will not burst at or below this pressure. Using a pressure monitoring device is recommended to prevent over pressurization. The size of stent expansion should be confirmed by radiography.

*Nominal Pressure: 10 atm (1013 kPa) for stent diameters up to 3.50 mm, 12 atm (1216 kPa) for 4.00 mm.

** Rated Burst Pressure: 18 atm (1824 kPa) for stent diameters up to 3.50 mm, 16 atm (1621 kPa) for 4.00 mm. **DO NOT EXCEED.**

14. REUSE PRECAUTION STATEMENT

Do not use if sterile barrier is damaged.
For single patient use only. Do not reuse, reprocess, or resterilize.

15. USE BY

The product storage period is valid for 12 months from the date of sterilization.

16. WARRANTY AND LIMITATIONS

Sino Medical Sciences Technology Inc. has committed itself to the manufacture of the HT Supreme™ stent system. Sino Medical Sciences Technology Inc. warrants that all products shall be free of defects in materials and manufacturing.

17. DECLARATION

Sino Medical Sciences Technology Inc. will not be liable for any accident and damage result from reuse of the product, incorrect size selection of the stent, misoperation or any other mistakes caused by human-factor (directly or indirectly).

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SINOMED B.V.