

RX

Tytrak™ PTCA Balloon Catheter

2.0 20

Tytrak™ 2.0 20

SN 201305199
2014 - 05 - 19
(YYYY - MM - DD)



(01)08717648118388(17)130131(10)102

Tytrak™ 2.0 20

SN 201305199
2014 - 05 - 19
(YYYY - MM - DD)

Model SNK-2.00-20-RX

PTCA Balloon Catheter / Cathéter à ballonnet pour PTCA / Catetere a palloncino per PTCA / Catéter con balón para ACTP / PTCA-Ballonkatheter / PTCA-ballonkatheter / Cateter Balão PTCA / PTCA氣球導管



Pressure (atm)	Diameter (mm)	
	(kPa)	
4	405	3.00
6	608	NP 3.09
8	811	3.24
10	1013	3.37
12	1216	3.51
14	1418	RBP 3.55
16	1621	3.59

Tytrak™ RX

2.0 mm x 20 mm

SN 201305199

SINO MED

Tytrak™ RX 2.0 mm x 20 mm



(01) 0 6948546 20421 3



(17)131209 (21) 211098003

SN 201305199



2014 - 05 - 19
(YYYY - MM - DD)

SINO MED

Tytrak™ RX

2.0 mm x 20 mm

SN 201305199

SINO MED

Tytrak™ RX 2.0 mm x 20 mm



(01) 0 6948546 20421 3



(17)131209 (21) 211098003

SN 201305199



2014 - 05 - 19
(YYYY - MM - DD)

SINO MED

Tytrak™ RX

2.0 mm x 20 mm

SN 201305199

SINO MED

EC REP

Alchimedics S.A.S.

Add: 9, avenue Paul Verlaine Grenoble Cedex 2 - France

Zip: 38030

Tel/Fax: 00 33 682 182 429



0.009.0063

CE 2460