



RX**HY SUPREME™** Drug Coated
Coronary Stent System**2.50 15****RX****HY SUPREME™** **2.50 15**REF 8.007.1812
2014 - 05 - 19
(YYYY - MM - DD)

0514 140519305023101

RX**HY SUPREME™** **2.50 15**

SN 305023101

2014 - 05 - 19
(YYYY - MM - DD)Drug Coated Coronary Stent System / Medikamentenbeschichtetes Koronarstentsystem / Sistema de stent coronario recubierto de fármaco /
Système d'endoprothèse coronaire à revêtement médicamenteux / Sistem Stent Koroner Berlapis Obat / Sistema stent coronario rivestito di
farmaco / Coronair stentsysteem met geneesmiddelcoating / Sistema de Stent Coronário Revestido por Fármaco / 藥物塗層冠狀動脈支架系統

	Pressure		Diameter
	(atm)	(kPa)	(mm)
8	811	2.38	
9	912	2.45	
10	1013	NP	2.50
11	1115	2.57	
12	1216	2.62	
13	1317	2.66	
14	1419	2.70	
15	1520	2.73	
16	1621	2.77	
17	1723	2.81	
18	1824	RBP	2.85

**HY SUPREME™** [RX] 2.50 mm x 15 mm

(01) 0 6948546 20514 1



(17)140519 (21) 305023101

SN 305023101 2014 - 05 - 19
(YYYY - MM - DD) **SINOMED****HY SUPREME™** [RX]

2.50 mm x 15 mm

SN 305023101

SINOMED**HY SUPREME™** [RX]

2.50 mm x 15 mm

SN 305023101

SINOMED**HY SUPREME™** [RX] 2.50 mm x 15 mm

(01) 0 6948546 20514 1



(17)140519 (21) 305023101

SN 305023101 2014 - 05 - 19
(YYYY - MM - DD) **SINOMED****HY SUPREME™** [RX]

2.50 mm x 15 mm

SN 305023101

SINOMED**EC REP**

SINOMED B.V.

Wilhelminakade 173, 3072 AP Rotterdam, The Netherlands

Tel: +31 10 307 6295

Fax: +31 10 307 6296



0.007.0115

CE 0344

339C

165*172mm